PROTOCOLS OF 406 MICROBIAL LABORATORY

Version2.0, Jan 2022

FREE for non-commercial use

Contributors

Jia-Min Kao, Shih-Han Wen, Chien-Yi Wu, Shu-Shuo Yeh,
Jing-Wen Michelle Wong, Emma Chen, Yan-Zhen Meng,
Hsuan-Tung Lin, Yen-Chih Lin

Contact

Shan-Hua Yang

TEL: +886-2-3366-2883

Email: shanhua@ntu.edu.tw

TABLE OF CONTENTS

1. SAMPLE PREPARATION

1.1.	Processing of coral samples.....	5
1.2.	Filtration of seawater samples.....	7
1.3.	Determination of Photosynthesis Activities.....	8
1.4.	Determination of Antioxidant Enzymes (SOD and CAT) Activities.....	9
1.5.	Cultivation	
1.5.1.	Formula of Marine Agar.....	19
1.5.2.	Formula of MMBV4 (Modified Marine Broth Version 4) Agar.....	21
1.5.3.	Formula of CHRO Agar (specify for <i>Acinetobacter</i>).....	23
1.5.4.	Formula of YM Agar (specify for yeast).....	24
1.6.	Preservation of bacteria.....	26
1.7.	Anaerobic culture apparatus and anaerobic culture.....	28

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1.3. Photosynthesis Determination (PAM 測定)

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1.4. Determination of Antioxidant Enzymes (SOD and CAT) Activities

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ROS (SOD, CAT) activities 實驗

B. Determine protein concentration in the sample

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Materials:

Kit content:

1. WST Solution
2. Enzyme Solution
3. Buffer Solution
1. Superoxide dismutase (SOD) Standard (Sigma, S7571-15KU)
2. 96-well microplate

Preparation of working solutions:

1. WST Working Solution: Dilute 1 ml of WST Solution with 19 ml of Buffer Solution. (Can storage at 4°C for 2 months, avoid light)
2. Enzyme working Solution: Centrifuge the Enzyme Solution tube for 5 seconds. Mix by pipeting, and dilute 15 μ l of Enzyme Solution with 2.5 ml of Dilution Buffer. (Can storage at 4°C for 3 weeks)
3. SOD Standard Solution (if necessary): Dilute SOD with Dilution Buffer to prepare SOD Standard Solution as follows:
 - a. 購得的 SOD standard 為粉末狀，總重僅 3 mg，以秤重方式無法配出精準濃度。故直接以 Dilution Buffer 加入藥瓶中回溶並稀釋成 1500 U/ml (stock)，每 150 μ l 分裝後冰入 -20°C 保存。取用時置於冰上退冰後使用 (只能單次使用，勿重複凍融)。

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




(4). Load 20 μ l Dilution buffer

(5). Load 20 μ l Enzyme working solution

*Avoid reaction time lag, use a multi-channel pipette

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(6) Incubate the plate at 37°C for 20 min

(7) Read the absorbance at 450 nm using a microplate reader

(8) Use sample absorbance to calculate the SOD activity (inhibition rate%)=

$$\frac{(Blank1-Blank3)-(Sample-Blank2)}{(Blank1-Blank3)} \times 100$$

Till this step, you should have: SOD standard concentration (A), SOD standard inhibition rate (B), Sample SOD inhibition rate (C)

(6) Calculate SOD calibration curve: reciprocal the SOD standard concentration (A) and its inhibition rate (B), draw a trendline base on these two parameters (1/A, 1/B) and get the trendline equation.

(7) Reciprocal the sample's SOD inhibition rate (1/C), then use the trendline equation to calculate sample's Reciprocal SOD concentration (1/D).

(8) Reciprocal back to get the sample's SOD concentration (D, U/ml), then divide it with the sample's protein concentration (Z, mg/ml) to get the sample's SOD concentration (U/mg)

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By Jing-Wen Michelle Wong

Materials:

1. Bacterial cultures
2. Plastic droppers
3. Glass slides

Solution


KOH (3%)	0.3 g
ddH ₂ O	10 ml

Methods:

1. Use inoculating loop to transfer generous amount of bacteria and spread on glass slide.
2. Apply 1 drop of KOH on the slide and stir carefully.
3. The solution of gram negative bacteria will be viscous and form a mucoid string within 30 sec.

Positive results: The solution with the bacteria (gram negative) will be viscous

Negative results: The solution with the bacteria (gram positive) will not be viscous

- A. Application of potassium hydroxide solution onto a microscopic slide.
- B. Bacteria are transferred and mixed with the KOH solution.
- C. Result with gram negative bacteria where the solution will be viscous and form a mucoid string.
- D. Result with gram positive bacteria where the solution will not be viscous.

2.2. Catalase Test

By Jing-Wen Michelle Wong

Materials:

1. Bacterial cultures
2. Plastic droppers
3. Glass slides

Solution


H ₂ O ₂ (1%)	1 ml
ddH ₂ O	9 ml

Methods:

1. Use inoculating loop to transfer generous amount of bacteria and spread on glass slide.
2. Add one drop of 3% H₂O₂ to the bacteria and observe the suspension (Be careful with the handling of H₂O₂ which is corrosive).

Positive test results: Gas formation (O₂) in the form bubbles shows that the bacterium has a catalase

Negative test results: No gas formation

- A. Apply one drop of 3% hydrogen peroxide on a microscopic slide.
- B. Transfer bacteria with a plastic loop to the H₂O₂ solution.
- C. Catalase positive bacteria produce gas (O₂) in form bubbles which shows that the bacterium has a catalase.
- D. Catalase negative bacteria does not have a catalase and therefore oxygen is not formed. Note the fragments of bacterial colonies, but there are no bubbles in the solutions.

2.3. Oxidase Test

By Jing-Wen Michelle Wong

Materials:

1. Bacterial cultures
2. Plastic droppers
3. Glass slides

Solution


TMPD (1.2%)	0.12 g
ddH ₂ O	10 ml

Methods:

1. Draw the filter paper into 4 sections, label with bacteria code.
2. Use inoculating loop to transfer generous amount of bacteria and spread on filter paper.
3. Apply 2 drops of oxidase reagent to the bacteria and observe it.

Positive test results: Dark blue-purple colour change within 10-30 sec.

Negative test results: No colour change /colour change after more than 30 sec.

- A. Two drops of oxidase reagent is applied onto a piece of filter paper.
B and C. Colony material is transferred by a plastic loop to the spot of oxidase on the filter paper.
C. The left spot includes oxidase positive bacteria and the spot to the right oxidase negative bacteria.
D. The final results can be observed after 30 s.
Note: Exposure to oxygen may eventually turn the test blue over time.

2.4. Lipase Test

By Jing-Wen Michelle Wong

Materials:

1. Bacterial cultures
2. Marine agar with Tween 80

Methods:

1. Draw the plate into 4 sections, label with bacteria code.
2. Use inoculating loop to transfer generous amount of bacteria and spread on plate.
3. Observe the plates after the bacterial colony grow.

Positive results: There are fatty droplets surround the colony.

Negative results: Without any droplets surround the colony.

Figure above shows positive results of lipase test.

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
File Index	Plate Num	Inc Hours	Plate Typ	Protocol	Strain Typ	Sample ID	Field 2	Field 3	Field 4	Field 5	Field 6	Field 7	Field 8	Field 9	Field 10
1	1	48	GEN III	A		IFA25-74-									
2	1	48	GEN III	A		IFA25-74-									
3	1	48	GEN III	A		IFA25-74-									

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


The screenshot shows the SoftMax Pro 7.1 software interface. The top menu bar includes tabs for Instrument, Temperature, Shake, Open / Close, Protocol Manager, Acquisition Settings, Read, Append, Reduction, Display, 3D, Zoom, Mask, Clone Plate, Template Editor, Copy, Import, Paste, Export, and Template Tools. The main window has tabs for Document, Comparison, and Workflow. The Navigation Tree on the left lists New Experiment, New Note, Delete Selection, BasicEndpoint, Intro, Expt1, Plate1, Standards, Unknowns, Unk_Dilution, Control, and StandardCurve. The central workspace displays a grid titled 'Plate1' with 12 columns labeled 1 through 12. The grid contains numerical values representing optical density measurements. A 'Reduction Settings' panel at the bottom shows 'Optical Density' and 'Wavelength Combination : !Lm1'. On the right, there are sections for 'Settings Information' (Endpoint Lm1 590, Shake Once, Calibrate On, Carriage Speed Normal, Column Priority) and 'Read Information' (SpectraMax ABSPlus ROM v2.0 20 Nov 05 2018, Start Read - 下午 02:24 2022/7/5). A message at the bottom right indicates a mean temperature of 28 °C.

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
Category

Wavelengths

- Plate Type
- Read Area
- PathCheck
- Shake
- Speed Read
- More Settings

Wavelength Settings

Number of Wavelengths Lm1 590 nm

You can specify up to six different wavelengths. Type the wavelengths in the fields.

[More Information](#)

Settings Information


- Endpoint
Lm1 590
- Read Area
Entire Plate
- Plate Type
96 Well Standard clrbtm
- More Settings
Shake Once
Calibrate On
Carriage Speed Normal
Column Priority

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
X

Read Mode: ABS

Read Type: Endpoint, Kinetic, Spectrum

Category: Wavelengths, Plate Type, Read Area, PathCheck, Shake, Speed Read, More Settings

Wavelength Settings: Number of Wavelengths: 1, Lm1: 750 nm

You can specify up to six different wavelengths. Type the wavelengths in the fields.

More Information

Settings Information:

- Endpoint: Lm1 750
- Read Area: Entire Plate
- Plate Type: 96 Well Standard clrbtm
- More Settings: Shake Once, Calibrate On, Carriage Speed Normal, Column Priority

OK Cancel

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3.1. Total DNA Isolations from Corals

3.1.1. Classical (CTAB) DNA Extraction

By Jia-Min Kao

改良自 “Current Protocols in Molecular Biology” 中“ Preparation of Genomic DNA from bacteria protocol”

Materials :

1. 1X TE buffer

Add 1 ml of 1M Tris-HCl (pH 7.5) in 98.8 ml MILIQ-water and then add 0.2 ml of 0.5M EDTA.

2. 10X TE buffer

Add 12.11 g of Tris-HCl (100mM) in 500 mL MILIQ-water and then add 3.722 g of EDTA (10mM). Adjust the pH up to 7.5 with NaOH(s), and then adjust volume to 1000mL.

3. 20 mg/ml proteinase K

4. Add 0.2 g of proteinase K in 10ml MILIQ-water.

5. RNase A

6. 10% SDS

Dissolve 20 g sodium dodecyl sulfate in MILIQ-water in a total volume of 200 ml with stirring.

7. 5M NaCl/6M NaCl

Dissolve 29.22 g (5M)/35.06 (6M) NaCl in 100ml MILIQ-water.

8. CTAB/NaCl solution

Dissolve 4.1 g NaCl in 80ml MILIQ-water and then add 10 g of CTAB with heating and stirring. Then adjust volume to 100 ml.

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4.1.3. Eukaryotic PCR Program (ITS1F & ITS4)

By Chien-Yi Wu

Fungus primer-PCR program (ITS1F & ITS4. Expected size:)

Create / Select a PCR program (machine ____ directory ____)

Stage	Step	Temp.	Dwell Time	Cycles
1	Step 01	94 °C	5 min	1
2	Step 01	94 °C	30 sec	
	Step 02	48 °C	30 sec	
	Step 03	72 °C	60 sec	
3	Step 01	72 °C	10 min	1
4	Step 01	10 °C	Hold	—

5.1. Nucleic Acid Purification

5.1.1. Gel Extraction by QIAEX II Gel Extraction Kit

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≤ 2 µg DNA	10 µl QIAEX II
2-10 µg DNA	30 µl QIAEX II
each additional 10 µg DNA	additional 30 µl QIAEX II

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By Emma Chen

Methods:

1. Add g sample to *PowerBead tube*, gently vortex to mix
 1. 0.3-0.5g for soft coral, ≤0.3 for hard coral
 2. For bacteria → a loop full
2. Add 60μl C1, invert several times to mix
 1. If C1 has precipitated, heat at 60°C until precipitate dissolves
3. Heat tubes at 70°C for 5min, then vortex for 5 sec → repeat 3x
 1. ***after this, set dry heat bath to 50°C*
 2. During this free time: can preemptively add C2, C3, C4 to Epps (see note at end of protocol)
4. Vortex tubes @ max speed (10) for 10min
5. Centrifuge @room temp: 10,000 x g, 30sec
 1. Small centrifuge takes time to get up to speed — if using this one, pick 10,000 x g, 1 min program (takes ~30 sec to get to 10,000)
 - i. In display: pick the correct program
 2. Large centrifuge: don't need to pick longer program
 - i. Can't pre-set programs, no need to check and adjust each time used. Also, only allows for times to be set in increments of 1 min, so set 30-sec protocols to 1min instead
6. Add 500μl of sample to 1.5ml Epp tube
 1. Take only clear part (not bottom/bubbles)
 2. Can toss PowerBead tubes after this
7. Add 250μl C2 to the tubes, vortex 5 sec (if not already added)
8. Put tubes in 4°C fridge for 5 min

1. Can leave them in there for ~30min-1hr if you wanna take a break — just write down how long you left them in fridge
9. Then centrifuge at room temp: 10,000 x g, 1min
10. Avoiding pellet, take 600µl supernatant → into 2nd 1.5ml tube
11. Add 200µl C3 (if not already added), vortex 5 sec, in 4°C fridge for 5min
 1. Same thing as pre-C2 addition*****
12. Centrifuge @room temp: 10,000 x g, 1 min
13. Avoiding pellet, take 750µl supernatant → into 2ml tube
 1. If not enough, just stop and DON'T take pellet
14. Add 1.2 ml C4 to 2ml Epp tube (if not already added), vortex 5 sec
 1. B/c it's over 1ml, just add 0.6µl twice
 2. **Shake C4 before use (good to shake other C solns too before use)
15. Load 675µl into a *SpinFilter*, centrifuge @room temp at 10,000 x g, 1min → DO THIS 3 TIMES
 1. After each time, empty the bottom level of spin filter
16. Add 500µl C5, centrifuge @room temp 10,000 x g, 30sec, then empty bottom level
17. Centrifuge @room temp 10,000 x g, 1min
18. Put SpinColumn into a *new 1.5ml tube*
 1. Avoid touching any remaining C5
19. Add 30µl ddH₂O right onto center of filter (right on top b/c there's just a tiny bit of it)
 1. Get ddH₂O from 4°C fridge (can be taken out earlier and put into 50°C to warm up)
20. Put them on 50°C heat dry bath for 30min
21. Centrifuge at room temp for 10,000 x g, 30 sec
22. Take out the SpinFilter & toss it, DNA should be in the bottom tube
23. Get 4 10µl Epp tubes per sample

1. Add 6.0μl of DNA sample into each 10μl Epp
 2. There's ~30μl per tube, so the last one will be more than 6μl (just put in the rest of it)
24. Store in -20°C or -80°C freezer

6.2. QIAEX II Agarose Gel Extraction (English Version)

By Emma Chen

Methods:

1. Take pics of all gels (GE samples & \oplus/\ominus control) → place back into plastic molds, bring to small room

In small room:

2. Slide gel out of plastic mold onto clear board, turn on light
 - a. Use orange shield to see gel bands (look from directly above to avoid cutting wrong part of gel)
3. **Cut each GE sample** (big band) out of gel
 - a. Use one disposable knife per GE sample → back into wrapper to dispose of it
 - b. Cut long sides first, then short sides, as close to the illuminated parts as possible
 - i. Make sure you cut deeply/fully so it comes out smoothly (w/o tearing)
4. Stab gel diagonally from the top long side to pick it up → place **into a 1.5ml Epp tube**

In hood/bench:

5. Add **700μl Q1** to each tube
 - a. Can take 1-2hr break after this step *only* if needed (samples → 4C fridge)*
6. Gently **vortex** the Q2 for ~30sec (due to its viscosity) → vortex @3.5
7. Add **10μl Q2** to each tube
 - a. Add this volume b/c we're working w/ $\leq 2\mu\text{g}$ DNA
8. Places tubes in floater divide in the **50C water bath** for **12 min** total
 - a. (official protocol = 10 min, JiaMin always does 12 min instead)

- b. Invert couple times every 1-2 min (in lieu of vortexing, as official protocol calls for)
 - c. After 12 min: make sure it's fully dissolved (not in bits) → dry off w/ paper towel
9. Centrifuge at 13,000 rpm / 2 min
- a. (official prot says 1min, JiaMin always does 2min for centrifuge steps)
 - b. Pipette out as much supernatant as possible (set pipette to ~700μl)
10. Add 500μl Q1, gently vortex
11. Centrifuge at 13,000 rpm / 2 min
- a. Pipette out as much supernatant as possible
12. Add 500μl PE Buffer, gently vortex
- a. PE Buffer comes as powder → 1st user must add ethanol (from kit, see volume on label), date it
13. Centrifuge at 13,000 rpm / 2 min
- a. Pipette out as much supernatant as possible
14. Add another 500ul PE Buffer
- a. ***DON'T VORTEX***this time!! → invert a couple times to wash pellet, BUT *ensure pellet's intact!!!*
15. Centrifuge at 13,000 rpm / 2 min
- a. Pipette out as much supernatant as possible (use blue then yellow pipettes to get it all out)
 - i. Will help air-drying proceed more quickly
16. Air-dry pellets: set tubes upside down on rack, place over airflow at the back of the hood for 25-30 min
- a. Place ddH₂O into 50C dry OR water bath to warm up for next step
 - b. ***Avoid over-drying!!!*** → set timer for 25 min, check → if not dry yet, wait another 5 min
 - i. When dry: pellets should be white (not translucent), w/ clearly defined edges

17. After dry: add 35µl ddH₂O
 - a. Can add 30µl or 35µl depending on which samples used (for QianYi' s → add 35µl)
18. Leave in hood for 30 min, gently vortex
 - a. Can leave at room temp b/c we're working w/ ≤ 4kb DNA (V6-V8 region = 424bp)
 - b. (official prot says 1min, JiaMin always does 30min) → BUT can leave ~1hr if need a break*
19. Centrifuge at 13,000 rpm / 2 min
 - a. DNA will be in the supernatant!!
20. End up w/ 20+µl supernatant, pipette it into two 10µl tubes, avoiding the pellet!!
 - a. One 10µl = 20µl, other 10µl = the rest of the volume
21. Measure DNA quality & concentration w/ NanoDrop

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Taiwan

<https://cran.csie.ntu.edu.tw/> ←

National Taiwan University, Taipei

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Download and Install R
Precompiled binary distributions of the base
one of these versions of R:

- [Download R for Linux](#)
- [Download R for \(Mac\) OS X](#)
- [Download R for Windows](#)

Subdirectories:

[base](#)

[contrib](#)

[old_contrib](#)

[Rtools](#)

3

[Download R 4.0.3 for Windows](#) (85 megabytes, 32/64 bit)
[Installation and other instructions](#)
[New features in this version](#)

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
Species	Treatment	Soil_type	Sampling_date
D001-16S	Diseased	Bulk_soil	T_0703
D002-16S	Diseased	Rhizome_soil	T_0703
D003-16S	Diseased	Bulk_soil	T_0703
D004-16S	Diseased	Rhizome_soil	T_0703
D005-16S	Diseased	Bulk_soil	T_0703
D006-16S	Diseased	Rhizome_soil	T_0703
D007-16S	Diseased	Bulk_soil	T_0710
D008-16S	Diseased	Rhizome_soil	T_0710
D009-16S	Diseased	Bulk_soil	T_0828
D010-16S	Diseased	Rhizome_soil	T_0828
D011-16S	Diseased	Bulk_soil	T_0904
D012-16S	Diseased	Rhizome_soil	T_0904
D013-16S	Diseased	Bulk_soil	T_0904
D014-16S	Diseased	Rhizome_soil	T_0904
D015-16S	Diseased	Bulk_soil	T_0911
D016-16S	Diseased	Rhizome_soil	T_0911
D017-16S	Diseased	Bulk_soil	T_0911

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```
#import taxonomy table
tax <- read.csv(file = "taxonomy.csv", row.names = 1, header = T, check.names = T)
tax <- as.matrix(tax)
TAX <- tax_table(tax)
```

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```
#filter data
#match the column and target, taxa in OTU table will decrease in every filter
filtered_physeq <- subset_taxa(physeq, Kingdom != "Unassigned")
filtered_physeq <- subset_taxa(filtered_physeq, Class != "Chloroplast")
filtered_physeq
```

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```
#merge same name in 1 in targrt taxonomy level(Phylum as example) for barplot, rest  
of the level will be "NA"  
filtered_physeq_abd_P <- tax_grom(filtered_physeq_abd, taxrank = rank_names  
(filtered_physeq_abd)[2])  
tax_table(filtered_physeq_abd_P)
```

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```
#####transfer sample count and rarefy####  
filtered_physeq_abd_P <- transform_sample_counts(filtered_physeq_abd_P, function(x)  
x / sum(x))  
  
filtered_physeq_abd_P_rarefy <- rarefy_even_depth(filtered_physeq_abd_P,  
sample.size = min(sample_sums(filtered_physeq_abd_P  
)),  
replace = F)
```

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The screenshot shows the measurement software interface with various windows and toolbars. A red box highlights the '測量' (Measurement) button in the top toolbar. Another red box highlights the '開啟計數器' (Open Counter) dialog box in the left sidebar, which lists three categories: '類1' (Category 1), '類2' (Category 2), and '類3' (Category 3). A third red box highlights the '另存新檔' (Save As New File) button in the top right corner of the main window. A fourth red box highlights the '導出到Excel表' (Export to Excel Table) button in the bottom right corner of the counter dialog. Step 1: Click on the image. Step 2: Enter the measurement page. Step 3: Open the counter (add category and check the category you want to calculate). Step 4: Click on the image to perform calculations. Step 5: Output results. Step 6: Save as new file.

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